

Patient Information Concerning NADAL® FOB


Store the sample collection tube in room temperature for usage (15 – 30 °C).

1. Gather a little bit of faeces using a stool collector or directly on toilet paper
Avoid diluting the faeces with water or urine.



2. Ensure that the sample collection tube is held upright, and open the BLUE cap.
Remove the enclosed sample collector.

3. Insert the end of the sample collector into the stool in three different places
consecutively. Faeces only in the grooves.



4. Place the sample collector with the stool sample in the sample collection
tube again, and close tight.

5. Note the date and time. Store the tube in a cool environment (2 – 8 °C)
where it can't be reached by direct light. Hand in the stool sample within
1 – 5 days.



If multiple samples have been ordered, these should be taken at separate times.